

## PLATING MEDIA FOR THE IDENTIFICATION OF *SALMONELLA*

The present invention relates to an isolating plating medium enabling the presumptive identification of *Salmonella*, and to methods of isolating bacteria of a particular type from other bacteria in a plating medium.

### BACKGROUND OF THE INVENTION

*Salmonella* is pathogenic for man, and the identification of *Salmonella* bacteria is a significant problem in medical bacteriology and in monitoring foodstuffs. It is of particular concern in the poultry and egg industries.

There have been many efforts to develop a selective plating medium for identifying *Salmonella* bacteria in samples containing mixed bacteria. "The Compendium of Methods for the Microbiological Examination of Foods" published by the American Public Health Association, 1992, describes some of these prior art processes at pages 382-383. In general, the prior art plating media for identifying *Salmonella* bacteria in a sample of mixed bacteria lacked specificity, were difficult to read, and were subject to false positives and/or false negatives. Further, many of the plating media of the prior art provided nutrients that are more advantageous for strains of bacteria other than *Salmonella*, hence providing a faulty indication of the presence of *Salmonella*. As a result there has been a continuing effort to develop a plating medium for identifying *Salmonella* in a mixed sample that overcomes or reduces the deficiencies of such media of the prior art.

The prior art includes Patent No. 4,279,995 dated July 21, 1981, to Woods et..al. entitled "Selective *Salmonella* Carbohydrate and Medium Constructed Therefrom"

which discloses a plating medium with 2-Deoxy-D-Ribose as a selective carbohydrate for *Salmonella* and a pH indicator dye to respond to carbohydrate metabolism. The media allows the growth of *Salmonella spp.*, *Arizona spp.* to the exclusion of other *Enterobacteriaceae*, but it also permits the growth of *Citrobacter freundii*.

Patent No. 5,098,832 dated March 24, 1992, and divisional Patent No. 5,194,374 to Rambach entitled "Isolating Medium for Identifying the *Salmonella* Bacterium" disclose a plating medium with 1,2-propanediol/silica gel that is metabolizable by *Salmonella* and a pH indicator to react to acidification of the medium. The patents also disclose the addition of a beta-galactosidase chromogenic substrate to increase the specificity of the medium.

Patent No. 5,434,056 dated July 18, 1995, to Monget et al., entitled "Method of Bacteriological Analysis, and Medium for the Detection of Bacteria of the *Salmonella* Genus" discloses a plating medium in which the acidic fermentation of sodium glucuronate monitored with a pH indicator, and a beta-D-galactopyranoside chromogenic substrate are used to facilitate selection of *Salmonella* colonies.

Patent No. 5,786,167 dated July 28, 1998, to Tuompo, et al., entitled "Method and Culture Medium for Identification of *Salmonella*" discloses a plating medium in which the acidic fermentation of melibiose, mannitol, and sorbitol is monitored with a pH indicator, and a beta-galactosidase responsive chromogenic substrate is used to facilitate selection of *Salmonella* colonies.

Patent No. 5,871,944 dated February 16, 1999, to Miller et.al., entitled "*Salmonella* Preferential Media" disclose a plating medium containing lactose and cellobiose as the carbohydrate source. In one embodiment, the production of H<sub>2</sub>S is used

to indicate *Salmonella*. In another embodiment, 5-bromo-4-chloro-3-indolyl-beta-D-galactopyranoside is incorporated in the medium and this chromogenic substrate is used to provide a color change to indicate bacteria other than *Salmonella*.

As indicated by the prior art set forth above, there continues to be a need for a plating medium for the differentiation of *Salmonella* that is more specific, easier to read, produces fewer false positives or false negatives, and is more sensitive to *Salmonella* bacteria.

#### SUMMARY OF THE INVENTION

The 1981 Patent No. 4,279,995 of Woods and Wilkinson, *supra*, discloses the process of selectively differentiating *Salmonella spp.*, *Arizona spp.* and some of the *Citrobacter freundii* from other members of the *Enterobacteriaceae* by using a medium containing 2-Deoxy-D-Ribose as a carbohydrate source and monitoring the metabolic acid of the medium. This patent further teaches the use of inhibitors to reduce or eliminate response of the medium to *Citrobacter freundii* and other non-target bacteria, but inhibitors also tend to have an adverse effect on the growth of target bacteria, particularly some strains of *Salmonella*. In accordance with the present invention, at least two chromogenic substrates are added to the differentiating medium to color colonies of non-target bacteria with essentially the same distinctive color, a color that contrasts with both the medium and the color produced by detection of the metabolic activity of *Salmonella* bacteria in the medium.

A plurality of chromogenic substrates is desirable to differentiate a sample containing a mixture of members of the *Enterobacteriaceae* because a positive reaction between the substrate and non-target bacteria is necessary, and chromogenic substrates

are generally selective. For example, *Citrobacter* does not react with beta-glucuronapryronaside but is reported to be 93% positive to beta-galactopyronoside. Hence, the present invention, in its preferred embodiment, utilizes beta-galactopyranoside substrates and effectively eliminates a response from *Citrobacter* which takes the color of the substrate.

Many of the bacteria that are found in mixed samples can be removed from the differentiation process by inhibitors without adversely effecting growth of *Salmonella* bacteria, particularly bacteria that are not member of the *Enterobacteriaceae*. Accordingly, a preferred embodiment of plating medium according to the present invention contains inhibitors.

#### DETAILED DESCRIPTION OF THE INVENTION

The selection of the carbohydrate, and the substrates determines the selectivity of a plating medium according to the present invention. A plating medium for differentiation of *Salmonella* bacteria, according to the present invention, can use most carbohydrates that produce metabolic reactions with *Salmonella*, but preferably the carbohydrate will not react with other bacteria, particularly other members of the family *Enterobacteriaceae*. 2-Deoxy-D-Ribose, xylose, mannitol, dulcitol, sorbitol, L-rhamnose and D-arabitol have been found to be suitable, and of this group of carbohydrates, 2-Deoxy-D-Ribose is preferred because of its strong positive reaction with *Salmonella* bacteria, including *Salmonella typhi*, and very few other bacteria of the *Enterobacteriaceae*.

The selection of the substrates depends upon the carbohydrate selected, since it will be the function of the substrates to change the color of colonies of non-target bacteria

to a common color that may be substantially ignored when assaying *Salmonella* colonies. The substrate may be either chromogenic or fluorogenic, and it is intended that the term substrate include both forms of substrates. The preferred substrates for differentiation of *Salmonella* are 5-bromo-4-chloro-3-indoxyl-beta-D-galactopyraniside, 5-bromo-6-chloro-3-indoxyl-beta-D-galactopyraniside, 3-indoxyl-beta-D-galactopyraniside, 6-chloro-3-indoxyl-beta-D-galactopyraniside, 4-nitrophenyl-beta-D-galactopyranoside, 2-nitrophenyl-beta-D-galactopyranoside, 5-iodo-3-indoxyl-beta-D-galactopyranoside, 4-methylumbelliferyl-beta-D-galactopyraniside and N-methylindoxyl-beta-D-galactopyranoside, but other substrates that do not react with *Salmonella* bacteria and do react with non-target bacteria may also be used. The substrates should be selected to produce water insoluble precipitate of approximately the same color from non-target bacteria.

It is also desirable to increase the production of precipitate from the chromogenic substrates by including an enhancer in the medium. The preferred embodiment of a plating medium according to the present invention includes the enhancer Isopropyl-beta-D-thiogalactopyranoside. Other suitable enhancers are 1-O-Methyl-beta-D-galactopyranoside, Ethyl-beta-D-thiogalactopyranoside, and Methyl-beta-D-thiogalactopyranoside

A pH indicator dye is required to produce *Salmonella* colonies of a contrasting color with the colonies of non-target bacteria. In the preferred embodiment of the invention, the non-target bacteria produce colonies of deep blue or purple, and a suitable pH indicator dye that will contrast with the deep blue may be red. Accordingly, the preferred embodiment incorporates a neutral red indicator dye. It is to be understood that

other dyes may be used, provided the dye contrasts with the color produced by the chromogenic substrates and with the color of the medium.

A number of inhibitors are desirable in the plating media of the present invention. One group of organisms that may readily be inhibited are gram-positive bacteria which are inhibited by bile salt and bile salt #3. Other inhibitors that may be used in the media of this invention are tellurite to retard the growth of *Escherichia*, sodium novobiocin to inhibit *Proteus* sp. and cefsulodin to inhibit *Pseudomonas/Aeromonas* sp.

The following table sets forth the ingredients of a plating medium that constitutes the preferred embodiment of the present invention.

TABLE 1

<u>Chemical</u>	<u>Grams/liter</u>
Yeast Extract	3.00
Proteose Peptone	10.00
Lab Lemco Powder	1.00
Sodium Chloride	5.00
L- Phenylalanine	3.50
Ferric Ammonium Citrate	0.50
Bile Salts #3	0.40
Bile Salts	0.20
5-bromo-4-chloro-3-indoxyl-beta-D-galactopyranoside	0.15
3-indoxyl-beta-D-galactopyranoside	0.15
Isopropyl-beta-D-thiogalactopyranoside	0.10
Neutral red	0.03
Agar	15.97
Deionized Water	950 ml/l.
<u>Supplements</u>	
2-Deoxy-D-Ribose	12.00
Sodium novobiocin	0.02
Cefsulodin	0.006

When preparing the medium, the ingredients, excepting the supplements, are mixed together in any order and thereafter boiled and cooled to form a basal medium. Thereafter, the supplements are added to the cooled, boiled basal medium just prior to completion of the plating medium.

Table 2 sets forth the results of tests made with a medium according to the preferred embodiment of the invention, as set forth in Table 1 above. The bacterial strains were incubated on the medium of Table 1 for 24 hours at 35 degrees Celsius, and Table 2 reports the results..

TABLE 2

Bacterial Strains	# of Strains	Colonial Morphology
<i>Salmonella spp.</i> (non <i>S. typhi</i> )	40	<ul style="list-style-type: none"> <li>• Raised colony; 1.5 - 4.0 mm in diameter,</li> <li>• Reddish pink color.</li> <li>• Some strains may have a pink precipitate surrounding the colony or a clear to light pink ring around the colony.</li> </ul>
<i>Salmonella typhi</i>	1	<ul style="list-style-type: none"> <li>• Raised colony; 1.0 - 1.5 mm in diameter.</li> <li>• Reddish pink color with no precipitate surrounding the colony</li> <li>• Clear to light pink ring around the colony.</li> </ul>
<i>Salmonella tennessee</i>	1	<ul style="list-style-type: none"> <li>• Raised colony; 2.0 - 3.0 mm in diameter.</li> <li>• Dark blue color with a pink precipitate and no clear ring around the colony.</li> </ul>
<i>Escherichia coli</i> 0157:H7	10	<ul style="list-style-type: none"> <li>• Domed to raised colony; 1.0 - 3.0 mm in Diameter,</li> <li>• Bluish green color with bluish green precipitate.</li> <li>• no clear ring around the colony.</li> </ul>
<i>Escherichia coli</i>	9	<p style="color: red;">e</p> <ul style="list-style-type: none"> <li>• Domed to raised colony; 2.0 - 4.0 mm in Diameter.</li> <li>• Bluish green color with bluish green precipitate.</li> <li>• No clear ring around the colony.</li> </ul>
<i>Escherichia coli</i>	4	<ul style="list-style-type: none"> <li>• Domed to raised colony; 2.0 - 3.0 mm in Diameter.</li> <li>• Dark blue in color with a pink precipitate.</li> <li>• No clear ring around the colony</li> </ul>
<i>Escherichia coli</i>	1	<ul style="list-style-type: none"> <li>• Raised colony; 2.0 mm in diameter..</li> <li>• Reddish pink color with a pink precipitate.</li> <li>• Thin clear ring around the colony.</li> </ul>
<i>Escherichia hermannii</i>	2	<ul style="list-style-type: none"> <li>• Domed colony; Pinpoint to 1.0 mm in Diameter.</li> <li>• Clear to greenish blue with no precipitate and ring around the colony.</li> </ul>
<i>Citrobacter diversus</i>	1	<ul style="list-style-type: none"> <li>• Domed to raised; Pinpoint to 1 mm in Diameter.</li> <li>• Greenish blue color with no precipitate.</li> </ul>
<i>Citrobacter freundii</i>	2	<ul style="list-style-type: none"> <li>• Domed colony; 2.0 - 3.0 mm in diameter</li> <li>• Dark blue to bluish green in color with no Precipitate and ring around the colony.</li> </ul>
<i>Serratia marcescens</i>	1	<ul style="list-style-type: none"> <li>• Domed colony; 2.0 mm in diameter.</li> <li>• Light green in color with no precipitate.</li> <li>• Thick clear ring around the colony.</li> </ul>

<i>Hafnia alvei</i>	4	<ul style="list-style-type: none"> <li>• Domed colony; Pinpoint to 1.0 mm in diameter.</li> <li>• Bluish green in color with bluish green precipitate; no ring around the colony</li> </ul>
<i>Enterobacter agglomerans</i> <i>Enterobacter cloacae</i> <i>Enterobacter aerogenes</i> <i>Enterobacter sakazakii</i>	4	<ul style="list-style-type: none"> <li>• Domed colony; 1.0 - 3.0 mm in diameter.</li> <li>• Bluish green in color with no Precipitate.</li> <li>• With or without a clear thin ring around the colony.</li> </ul>
<i>Klebsiella ozaenae</i>	1	<ul style="list-style-type: none"> <li>• Domed colony; 1.0 - 2.0 mm in diameter.</li> <li>• Clear to tan in color with no precipitate.</li> </ul>
<i>Klebsiella pneumoniae</i>	2	<ul style="list-style-type: none"> <li>• Domed colony; 2.0 - 3.0 mm in diameter.</li> <li>• Bluish green in color with bluish precipitate arid no clear ring around the colony.</li> </ul>
<i>Morganella morganii</i>	1	<ul style="list-style-type: none"> <li>• Domed colony; 1.0 - 2.0 mm in diameter.</li> <li>• Clear to cream color with a brownish precipitate in the medium. No clear ring around the colony.</li> </ul>
<i>Providencia retigeri</i> <i>Providencia alealifaeiens</i> <i>Providencia stuartii</i>	3	<ul style="list-style-type: none"> <li>• Domed colony; Pinpoint to 2.0 mm iii diameter.</li> <li>• Clear to tan color with brownish Precipitate in the medium. No clear ring around the colony.</li> </ul>
<i>Acineobacter calcoueticus</i> <i>Proteus mirabilis</i> <i>Pseudomonas aeruginosa</i> <i>Yersinia enterocolitica</i> <i>Pseudomonas pickettii</i> <i>Aeromonas hydrophila</i>	8	<ul style="list-style-type: none"> <li>• No growth for all strains tested.</li> </ul>

From Table 2, it is clear that the medium of the present invention is highly selective for *Salmonella spp.* and *Salmonella typhi* bacteria, and that samples containing mixed bacteria on plates of that medium that have been properly incubated are readily assayed because the target colonies are uniquely colored. The use of a carbohydrate and a pH indicator dye to color colonies of the target bacteria and multiple chromogenic substrates to color colonies of non-target bacteria has provided improved selectivity and facilitated assaying. *Salmonella* bacteria are positive with respect to the selected carbohydrate, and negative with respect to the substrates, thus forming colonies with the

color of the dye. Non-target bacteria are positive with respect to one or both substrates. If the non-target bacteria are negative with respect to the carbohydrate, colonies of these bacteria assume the color of the substrate. If the non-target bacteria is also positive with respect to the carbohydrate, colonies of these bacteria will assume a color that is a blend of the color of the dye and the active substrate.

While the plates of the preferred embodiment of this invention are designed for the identification of *Salmonella* bacteria, those skilled in the art will be able to adapt this invention readily for the identification of other microorganisms. It is therefore intended that the scope of this invention be not limited by the specification, but rather only by the appended claims.

098855204-0620014